

Clavicle Plating System --- Set Configuration

Bone Plates, CP titanium			
Part Number	Length (mm)	Holes	Qty.
601-06	72	6	2
601-07	84	7	2
601-08	97	8	2
601-09	110	9	2
601-10	123	10	2
601-11	136	11	2



Case	
Part Number	Qty.
900-01	1

4mm Bone Screws, titanium alloy		
Part Number	Length (mm)	Qty.
402-12	12	10
402-14	14	20
402-16	16	20
402-18	18	10
402-20	20	10
402-22	22	10
402-24	24	5



2.5mm Fragment Screws, titanium alloy		
Part Number	Length (mm)	Qty.
251-10	10	5
251-12	12	5
251-14	14	5
251-16	16	5
251-18	18	5
251-20	20	5
251-22	22	5
251-24	24	5



Instruments		
Part Number	Description	Qty.
801-01	Clavicle Retractor	1
801-02	Hohmann Retractor	2
802-01	Bone Reduction Forceps	2
802-02	Stagbeetle Forceps	1
803-01	Bending Iron	2
804-01	Depth Gauge	1
805-01	Drill Guide	1
806-01	Palm Handle	1
806-02	Jeweler's Handle	1
807-01	1.5mm Hex Screwdriver Shaft	1
807-02	2.5mm Hex Screwdriver Shaft, Short	1
807-03	2.5mm Hex Screwdriver Shaft	1
808-01	1.9mm Drill Bit	2
808-02	2.4mm Drill Bit	2
809-01	2.5mm Fragment Screw Tap	1
809-02	4mm Bone Screw Tap	1
810-01	16° Countersink	1
810-02	Screw Retaining Sleeve	1



Medical, LLC

52 Highland Avenue, Suite B
Bethlehem, PA 18017
T. 610.954.7901
F. 610.954.7904
www.mbdmedical.com

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.